




LiteGait® LG50KXFD46



Max. Patient Weight	50 kg	Base Options	Standard Base	64	64/69/76	
Max. Patient Height Over Ground		Base Dimensions		64	69	76
Over Ground	160	Treadmill Width		62	67	74
Max. Patient Height w/Treadmill		Door Width		75	80	87
Over Treadmill	145	Max. Base Footprint		84 x 82		
Yoke Adjustment		Clearance		39		
Powered	45	Leg Length		82		
Manual	70	Height (Floor to Top of Leg)		20		
Overall Unit Height	178	Caster Diameter		12.7		
		Shipping Dimensions		76 x 46 x 119 cm		
		Shipping Weight		44 kg	45 kg	46 kg
		Power Requirements:		110V-240V 50/60Hz		

LiteGait® Includes			
Integrated Freedom	INCLUDED	Included Harnesses	2 Harness
BiSym	NOT AVAILABLE	Diaper	
Gaiter Stool	OPTIONAL	Pediatric	
3D Handlebar	NOT AVAILABLE		
Parallel Grab Bar	INCLUDED	Q-Straps	OPTIONAL

Customization Available - Contact Mobility Research Sales Dept. for additional information
Created: 12/12/2018

All dimensions subject to change
Dimensions listed in cm/kg
Treadmill step up height: 16 cm