


PRODUCT SHEET

LiteGait® LG900



Max. Patient Weight	225 kg
Max. Patient Height Over Ground	
Over Ground	230
Max. Patient Height w/Treadmill	
Over Treadmill	215
Yoke Adjustment	
Powered	90
Manual	N/A
Overall Unit Height	255

Base Options	Standard Base	81	76/81/86
Base Dimensions			76 81 86
	Treadmill Width		74 79 84
	Door Width		90 95 100
	Max. Base Footprint		97 x 122
	Clearance		52
	Leg Length		122
	Height (Floor to Top of Leg)		20
	Caster Diameter		15.2
Shipping Dimensions			107 x 71 x 140 cm
Shipping Weight			125 kg 130 kg 135 kg
Power Requirements: 110V-240V 50/60Hz			

LiteGait® Includes

Integrated Freedom	INCLUDED	Included Harnesses	2 Harness
BiSym	WIRELESS	Small Adult	
Gaiter Stool	OPTIONAL	Standard Adult	
		Ext. Panel w/43 cm Groin Piece	
3D Handlebar	INCLUDED		
Parallel Grab Bar	NOT AVAILABLE	Q-Straps	INCLUDED

Customization Available - Contact Mobility Research Sales Dept.
for additional information
Created: 12/12/2018

All dimensions subject to change
Dimensions listed in cm/kg
Treadmill step up height: 16 cm